


4. How to activate SOUL 2.0 Product?

Ans: There are 2 kinds of SOUL 2.0 keys: One is that Product Key a.k.a ID Key written on CD-Cover beneath Institute Name. This key is permanently provided and is used when SOUL 2.0 Server is installed

Second key is Handshake key that is generated after SOUL 2.0 server is installed and is required to provide to INFLIBNET Centre.

Third is Activation key that will be provided to the Institute by INFLIBNET Centre in return of the Handshake Key. Both Handshake Key and Activation key are required to activate SOUL 2.0 Server

- Once SOUL 2.0 Server is activated it will show the below given Screen



- Click on Activate Your Product which will generate the Handshake Key that has to be given through Email to SOULHELP0007@GMAIL.COM or by Fax 079-26300990
- Enter Activation Key in the Activation Key field and Click on Activate